

# SediGraph™ III 5120

Particle Size Analyzer

 micromeritics  
SediGraph III  
Particle Size Analyzer

# New Improvements to a Classic Technique

## SediGraph III 5120 Particle Size Analyzer

### Proven Technique and Reliability

For over three decades, the Micromeritics SediGraph has remained the standard instrument for particle size analysis in many laboratories throughout the world. Whether in a rugged production environment or a controlled laboratory setting, the SediGraph continues to produce accurate results with superior reliability. Particle size distribution is measured using the sedimentation method. Particle mass is measured directly via X-ray absorption. By measuring the rate at which particles fall under gravity through a liquid having known properties as described by Stokes' law, the SediGraph determines the equivalent spherical diameter of particles ranging from 300 to 0.1 micrometers.

The new generation SediGraph III 5120 combines this proven technique with new technology to deliver reproducible and highly accurate particle size information, completing most analyses in minutes.

### Intelligent Design Features

The SediGraph III 5120 offers advanced instrumentation features that ensure measurements are repeatable and easy to perform. New features make it easier to operate and maintain the instrument, and the results can be reliably reproduced by SediGraphs in other locations. Design improvements include:

- Utilization of a simplified pumping system ensures **fast and easy maintenance**
- **Reduced noise level** for a quieter working environment
- **A maintenance reminder**, based on the number of analyses performed, alerts you when it is time for routine maintenance
- **Computer-controlled mixing chamber temperature** improves repeatability and reproducibility
- **A highly versatile and interactive reporting system** provides a wide range of custom data presentation options and now includes particle settling velocity and grain size in Phi units



## A Wide Variety of SediGraph III 5120 Benefits

- **Complete particle accountability** assures that all of the introduced sample is accounted for, including any fraction above 300  $\mu\text{m}$  and below 0.1  $\mu\text{m}$
- **Capability to merge data** with that from other particle sizing methods, thus extending the range of reported data to 125,000  $\mu\text{m}$  (125 mm), excellent for geological applications
- **Scanning the sedimentation cell** from bottom to top allows accurate inventory of fast-settling particles while minimizing the time required to resolve the separation of fine particles
- **Fully automatic operation** increases sample throughput and reduces operator involvement in addition to reducing the opportunity for human error
- **Temperature-controlled analyses** assure that liquid properties remain constant throughout the analysis so you can be confident of accurate results
- **Multiple analysis speeds** allow you to choose the desired combination of speed and resolution that meets your needs
- **Real-Time display** allows you to monitor the cumulative mass plot of the current analysis and to make immediate procedural changes if needed
- **Statistical process control (SPC) reports** track the performance of your processes allowing immediate response to fluctuations
- **Plot overlays** provide a visual comparison of analysis results from one or more analyses; a reference or baseline analysis, for example, or a superposition of two different types of plots of the same analysis data
- **Data comparison plots** provide graphical displays of the mathematical difference between two data sets (difference from reference plot) or the extent of a data point value above or below a tolerance boundary (out of specification plot)
- **Multiple analyzer control** allows two SediGraph III's to be operated simultaneously from a single computer, conserving valuable lab space and making data storage convenient

## Applications

**Ceramics:** The size range of particles and the distribution of mass in each size class strongly affect the ability to sinter a ceramic powder and its forming properties as well as the pore size distribution in the finished product. Particle size distribution information helps determine curing and bonding procedures, control pore structure, ensure adequate green body strength, and produce a final product of desired strength, texture, appearance, and density.

**Metal Powders:** By controlling particle size, very specific pore characteristics can be designed into a product. Porosity characteristics often are the key to product performance. Similar to ceramics, the particle size distribution is critical to green body and final product strength and density.

**Geological/Soil Science:** Grain size affects the moisture-holding capacity of soil, drainage rate, and the soil's ability to hold nutrients. Grain size is directly related to transport of sediment.

**Cosmetics:** The appearance, application, and packaging of cosmetics are influenced by the particle size distribution of base powders, such as talc, and the pigments used for coloring.

### Pigments:

Particle size alone can affect the tinting strength of a color. As tinting strength goes up, the quantity of pigment needed to produce required color intensity goes down. The particle size affects the hiding power of the paints. Also the particle size distribution influences gloss, texture, color saturation and brightness.

**Catalysts:** Particle size affects the catalytic activity of a metal for structure-sensitive catalytic reactions.

**Construction Materials:** Particle size of cement affects setting time and strength characteristics of the finished concrete and cement.

**Minerals and Inorganic Chemicals:** Reactivity of materials is dependent upon exposed surface area and thus particle size distribution.


**Abrasives:** A properly balanced size distribution of abrasive grains and powders is a fundamental consideration whether the material is to be used in slurries, dry blasting, or bonded abrasive tools. Uniform particle size assures precise flow rates through blast machines and is a critical determination in media management when recycling the abrasive material.



# An Enduring and Proven Particle Size Method

## The Endurance of the SediGraph Method


The SediGraph has been employed in a wide variety of industrial applications since its embodiment into a commercial instrument in 1967. To confirm its wide use in various applications throughout the world, one needs only to type in 'sedigraph' as the search key for any Internet search engine. The instrument has undergone many improvements in speed, sample handling, and data reduction and reporting since its introduction. However, the fundamental analytical technique continues to be based on two well-established and well-understood physical phenomena—sedimentation and photon absorption. Stokes' law is applied to determine particle size by measurement of the terminal settling velocities of sample particles of various sizes. Relative mass concentration for each size class is determined by applying the Beer-Lambert-Bouguer law to the measured absorption of a low-power X-ray beam projected through the fraction of sample remaining in suspension. The elegant simplicity of the Stokes and Beer-Lambert-Bouguer laws means that interpretation of raw data is straight-forward; the analyst easily can understand the relationship between the basic measurements and the reported size distribution. All experimental parameters are easily determined, data reduction is uncomplicated and fast, and there is no requirement to 'bias' the data reduction software toward a particular distribution modality.



## Stokes' Law

Stokes' law simply states that the terminal settling velocity of a spherical particle in a fluid medium is proportional to the square of the diameter of the particle. Stokes' law applies rigorously providing that a certain relationship between these variables, the particle diameter and the settling velocity, is not violated. Application: if a collection of particles of various diameters is uniformly dispersed in a liquid of density less than that of the particles and then allowed to settle under gravity, it accurately can be predicted when all particles below any given size will have fallen below any given level. The size distribution of the particles can be extracted from this information. For a detailed description of the SediGraph technique refer to ISO 13317-3:2001 Determination of particle size distribution by gravitational liquid sedimentation methods – Part 3: X-ray gravitational technique.






## Quality Measurement Results

Critical decisions or conclusions may be based on size measurements, therefore your confidence in the quality of those measurements should be well-founded. Data *accuracy* is particularly important in research applications. Accuracy describes how closely the measurement agrees with the accepted or 'true' value. *Repeatability* and *reproducibility* are of equal importance in quality control applications. Repeatability describes the ability to produce the same results over a series of measurements of the same sample. Reproducibility of measurements usually pertains to different instruments in different facilities, an important instrument quality in assuring that each facility will produce equivalent products.

An easy-to-use instrument with the capability to produce precise, repeatable, and reproducible data reliably and automatically at high throughput rates epitomizes the ideal measuring tool. The SediGraph is renowned for repeatability and reproducibility and is used as a quality control instrument in many industries that require consistent product to be produced in every factory location. Likewise, the accuracy of the SediGraph also has earned it a high level of respect. The SediGraph directly measures the X-ray absorption mass of the sample and reports mass distribution – no modeling is required! Many particle size reference materials used to test sizing instruments are qualified by the manufacturer using a SediGraph.

## Sedimentation, X-ray Absorption

The SediGraph uses a narrow collimated beam of X-rays to measure directly the particle concentration in the liquid medium. This is done by first measuring the intensity of a baseline or reference X-ray beam which is projected through the cell windows and through the liquid medium prior to the introduction of the sample. A homogeneously dispersed mixture of solid sample and liquid is next pumped through the cell. The attenuated X-ray beam is measured to establish a value for full scale attenuation. Agitation of the mixture is ceased and the dispersion is allowed to settle while X-ray intensity is monitored. During the sedimentation process, the largest particles fall below the measuring level, and progressively finer and finer particles do so until only the finest remain near the top of the measuring cell.



# Extremely Versatile Data Presentation and Reporting System

The SediGraph III 5120 is equipped with a versatile, easy-to-use user interface that provides all of the convenient features you expect from a Windows®-based program. These features include point-and-click menus, customizable reports with your laboratory logo graphic, editable graphs, cut-and-paste graphics and tables, data export features, and more. Custom protocols help plan, launch, and control the analysis and assure that subsequent analyses are all performed in the same manner, regardless of the skill of the operator. You can collect, organize, archive and reduce raw data, and store standardized sample information and analysis conditions for easy access during later applications. Finished reports may be generated to screen, paper, or transferred in a variety of formats to storage devices.

Report by Size Table				
High Diameter (µm)	Low Diameter (µm)	Low Diameter (φ)	Cumulative Mass Finer (Percent)	Settling Velocity (cm/s)
250.0	200.0	2.322	51.2	4.79129
200.0	150.0	2.737	43.3	2.69510
150.0	100.0	3.322	38.4	1.19782
100.0	80.00	3.644	33.2	0.76661
80.00	60.00	4.059	30.2	0.43122
60.00	50.00	4.322	27.9	0.29946
50.00	40.00	4.644	24.5	0.19165
40.00	30.00	5.059	19.7	0.10780
30.00	25.00	5.322	15.7	0.07486
25.00	20.00	5.644	7.3	0.04791
20.00	15.00	6.059	4.5	0.02695
15.00	10.00	6.644	3.5	0.01198
10.00	8.000	6.966	2.9	0.00767
8.000	6.000	7.381	2.2	0.00431
6.000	5.000	7.644	1.9	0.00299
5.000	4.000	7.966	1.5	0.00192
4.000	3.000	8.381	0.9	0.00108
3.000	2.000	8.966	0.0	0.00048


## Data Reporting

Detailed analysis data for particles ranging from 300 to 0.1 µm are provided automatically by the SediGraph III. Data collected from other particle size analyses ranging from 125,000 to 300 µm can be combined with SediGraph data, enabling effective reporting for particles ranging from 125,000 to 0.1 µm. That fraction finer than 0.1 µm is also indicated.

In addition to tabular data, different graphical analysis plot types are available including:


- Cumulative Mass, Area, and Number
- Settling Velocity Distribution
- Process Control Charts
- Log Probability
- Baseline/Full Scale References
- Frequency Distribution
- Difference From Reference
- Out of Specification
- Rosin-Rammler
- Regression Analyses

*Tabular report and cumulative finer mass versus phi plot for soil sample*



Plots can be overlaid for comparing the results from different samples or for comparing different plot types from the same sample. This allows you to compare analysis results to a standard. Plots can be rescaled to give you the ability to scrutinize closely your graphical data.

A new column has been added to the tables and to the x-axis selection in the graphs that reports size in Phi units, where  $\Phi = -\log_2$  (particle diameter in mm). Also, a column selection for settling velocity (cm/s) is available for tables. The x-axis of plots can be scaled in particle size or settling velocity. The SediGraph is the only automated particle sizing instrument that directly measures settling velocity.



## SPC Reporting and Regression Analyses

Statistical Process Control (SPC) reporting provides an easy method for continuously monitoring production processes and reducing response times to deviations from the standard. Regression analyses allow you to determine the relationship between a control parameter, for example, and a measured characteristic of the sample. Select from 26 axis variables including:

- Mean
- Mode
- Median
- Size at (percentile)
- Standard Deviation
- Coefficient of Variation
- - N  $\sigma$  Size
- + N  $\sigma$  Size
- Skewness
- Kurtosis
- Specific Surface Area
- Cumulative Percent at Size
- Percent Out of Specification
- Full Scale Scan Pump Speed
- MasterTech Stirrer Speed
- MasterTech Stirrer Time
- MasterTech Ultrasonic Probe Time
- Particle Density
- Liquid Viscosity
- Liquid Density
- Three User-defined (External) Parameters

## Accessories

### The MasterTech™ 052 Autosampler

The MasterTech 052 Autosampler provides assurance that samples are prepared and analyzed exactly the same way, every time. The MasterTech is designed to increase throughput, repeatability, and reproducibility while reducing operator involvement. Up to 18 samples can be queued to run sequentially and completely unattended, including automatic stirring or sonication prior to transfer to the analysis system. The SediGraph III's operating software controls the MasterTech, and information about dispersion is stored in the sample file for future reference.

The MasterTech features a powerful ultrasonic probe for sample redispersion. Power to the probe tip is adjustable and the driving


circuit is self-tuning for maintaining efficient and consistent sonic energy levels. A front-panel digital readout lets you know when the desired power is reached, and that same power is applied each time the method is repeated. The SediGraph III 5120 is designed to accommodate the MasterTech on top of the analytical unit, thus conserving valuable bench space.



### SediSperse®

Homogeneous sample dispersion is crucial for obtaining accurate and reproducible results in any particle size study. Micromeritics produces a complete line of SediSperse dispersing liquids specifically designed to maximize particle suspension and eliminate agglomeration of particles. SediSperse liquids are available in nine variations of aqueous and organic formulations and are effective in dispersing most materials.

*Cumulative finer mass versus diameter plot and SPC control chart for fine calcium carbonate sample*



To request a quote or additional product information, visit Micromeritics' web site at [www.micromeritics.com](http://www.micromeritics.com), contact your local Micromeritics sales representative, or our Customer Service Department at (770) 662-3636.



*The Science and Technology of Small Particles™*

**Micromeritics Instrument Corporation**  
4356 Communications Drive  
Norcross, GA 30093  
USA

Telephones:

U.S.Sales (770) 662-3633  
International Sales (770) 662-3660  
Fax (770) 662-3696

**Micromeritics China - Beijing Office**  
Room 1202. No. 1 Building  
Shi-Hao Ming-Di (Ten-Mansion)  
No. 81 Zi Zhu Yuan Rd.  
Hai Dian District  
Beijing, P.R. CHINA  
Code: 100089

Telephone: (86)-10-68489371  
Fax: (86)-10-68489372

**Micromeritics France S.A.**

Parc Alata  
Rue Antoine Laurent Lavoisier  
F-60550 Verneuil en Halatte  
FRANCE

Telephone (+33) (0)3 44 64 60 80  
Fax (+33) (0)3 44 64 60 89

**Micromeritics GmbH**  
Avantis Science Park  
Rutherford 108  
D-52072 Aachen

Telephone (+49) (0) 241 189 446 0  
Fax (+49) (0) 241 189 446 11

**Micromeritics Ltd.**  
Unit 2, Chestnut House  
178-182 High Street North  
Dunstable, Bedfordshire LU6 1AT  
ENGLAND

Telephone (+44) (0)1582-475248  
Fax (+44) (0)1582-475252

**Micromeritics N.V./S.A.**  
Eugene Plasky laan 140B  
1030 Brussels  
BELGIUM

Telephone (+32) (0)2-743-39-74  
Fax (+32) (0)2-743-39-79

**Micromeritics SRL**  
Via W.Tobagi n. 26/7  
20068 Peschiera Borromeo, Milano  
ITALY

Telephone (+39) (0)255302833  
Fax (+39) (0)2553 02843

**Micromeritics Japan, G.K.**  
5F Tokatsu Techno Plaza  
501, 5-4-6 Kashiwanoha  
Kashiwa, Chiba 277-0882  
Japan

Telephone (81)-0-4-7128-5051  
Fax (81)-0-4-7128-5054