

Use of Minimum Signal Fraction Setting in Saturn DigiSizer Sample Information Files

The light scattered by particles is only a portion of the total laser light scattered during analysis with the Saturn DigiSizer. The optical components of the analyzer and the dispersing liquid also scatter light; this is known as *background scattering*. The light scattered by the particles compared to that of the background is called the *scattering signal fraction*. In order to limit the effects of scattering noise, a threshold known as the **Minimum signal fraction** can be specified in the Saturn DigiSizer sample information file. Any collected data that do not exceed the background by this specified intensity percentage are not used in particle size calculations.

The minimum signal fraction default value is 16%. This means that scattering data collected which have intensities equal to or greater than 16% above the background level at that angle are used in the particle size calculations. Data points where the scattering intensity is less than 16% above the background are not used in calculations, although data are collected at all angles where the intensity is 3% above the background.

If you need to increase the signal fraction at particular angles, increase the concentration of the sample. This results in a linear increase in

signal fraction and, therefore, data at more angles (generally smaller angles) are used. This also allows detection of larger particles which may have been missed in the original experiment. There is a limit, however, in using higher concentrations. If small particles are present, multiple scattering may result. In this case, what do you do in order to size the larger particles in the presence of the smaller ones in a sample with a broad particle size distribution? You change the minimum signal fraction to allow more data to be used.

In Figure 1 the original default minimum signal fraction of 16% is used to calculate the size of an alumina sample analyzed with an obscuration of 10.9%. The large end of the distribution is artificially truncated because the signal from the larger particles is less than 16% above background. For this sample, the concentration could not be increased to raise the signal for the larger

particles since doing so would result in multiple scattering by the smaller particles. This can be seen in Figure 2 when the sample is analyzed with an obscuration of 29.0%. In Figure 3 the minimum signal fraction is reduced to 5% for the original sample. Notice that the full distribution is calculated since the largest size is the same as when the concentration was increased (Figure 2).

So, should you always use a lower signal fraction rather than a higher concentration for samples that contain only large particles? No; use a higher concentration when possible. This improves the statistics, and thus the reliability, of the analysis since a larger number of total particles is analyzed. It also reduces the effects that noise can have on the analysis.

—mi






Figure 1. Alumina sample analyzed at 10.9% obscuration and 16% minimum signal fraction.




Figure 2. Alumina sample analyzed at 29.0% obscuration and 16% minimum signal fraction.




Figure 3. Alumina sample analyzed at 10.9% obscuration and 5% minimum signal fraction.

Micromeritics' AquaPrep™ Helps Eliminate Air Bubbles

The use of water as a dispersing liquid for particle size analysis can often result in misleading or even incorrect data due to the presence of air bubbles. If the air bubbles are not removed, the particle size analyzer may detect them and report them incorrectly in the measured distribution. To solve this problem, Micromeritics has developed the AquaPrep.

The AquaPrep operates by recirculating water through a hydrophobic capsule consisting of many thin-walled capillaries. A vacuum pump in turn provides suction on the outside of the capillaries. This results in a diffusion of dissolved air from the water through the capillary walls and into the vacuum pump. The air removed from the water is exhausted through a small tube at the front of the instrument. The AquaPrep is convenient, small in size, and simple to use. For more information on the web, see www.micromeritics.com/pdf/products/AquaPrep.pdf