








- The standard unit can be adapted for different catalytic tests through different configurations and options.
- Real time results with very high level of accuracy and reproducibility.
- Programming of series of experiment with data acquisition and graphics capabilities.
- Up to 1100 °C, depending on the material of reactors.
- Up to 200 bar as optional.


STANDARD SPECIFICATIONS

- 3 x high precision mass flow controllers with digital communications.
- Preinstallation of liquid feed line.
- 1 SS316 tubular reactor, 9.1 mm ID with easy loading of up to 3.3 cc catalyst. Thermocouple placed directly in catalyst bed.
- Maximum working pressure up to 100 ± 0.1 bar based on micrometric servo-controlled valve design. PED 2014/68/UE, high-pressure certification.
- Maximum working temperature up to 1100 ± 1 °C (depending on reactor material, SS316 as standard. 9,1 mm ID x 300 mm length, 20 μm porous plate.
- All layout inside hot box made with hot air convector. Maximum temperature up to 200 ± 1 °C.
- Automatic six-port reactor bypass valve.
- High pressure Liquid/Gas separator (L1) with capacitive level sensor and very low dead volume (less than 1 cc), allowing real time results without accumulation. Level control valve (LCV) based on micrometric servo-controlled valve.
- Independent safety levels separate from PC. User-defined functions for alarms.
- PC with user-friendly software with real time supervision and program recipes.



TOUCH SCREEN





4

PROCESS® SOFTWARE AND EXPERIMENT VIEWER

- Process® is a user-friendly real-time SCADA software based on LabVIEW and digital communications.
- Session table for automatic devices changes. Easy programming of conditional lines for more efficient experiments.
- Monitoring and data acquisition for process and control values.
- Graphs capabilities.
- All parameters are saved in a log file for later graphical representation and data treatment. Easy export to spreadsheets.






REACTOR MATERIAL AND SIZES

	Inside x Outside diameter: 9.1x14.3 mm (Std), 13.1x19.1 mm*, 17.5x25.4 mm*					
	Max. temperature (°C)	Max. pressure (BAR)	Max. P (BAR) @ 500 °C	Max. P (BAR) @ 650 °C	Max. P (BAR) @ 790 °C	
SS316 (Std)	800	690	214	104		
Hastelloy X*	1200	617	525	308	100	
Hastelloy C276*	1093	418	364	155		
Inconel 600*	1212	356	243	36		
Inconel 625*	980	827	728	414		
SS310*	1100	441	262	241		

23.8 mm ID reactor available

* Optional Std: standard




01 MICROACTIVITY EFFI

OPTIONAL CONFIGURATION

GTL configuration


High pressure L/L/G separator with very low dead volume (less than 1 cc each phase) for separation and collection of 2 non miscible liquid products.





2 reactors


1_Duo configuration

- 2 fixed bed reactors.
- Possibility to switch between reactor 1 or reactor 2 and serial connection.
- Optionally a 2nd PCV can make the reactor 1 work at higher pressure than reactor 2.
- Optionally an additional MFC fed to reactor 2 to modify the composition coming from reactor 1.










2_Twin configuration


- ▶ 2 fixed bed reactors, one for reaction for pressures up 100 bar, one for regeneration at atmospheric pressure.
- 1 x MFC for regeneration.
- Possibility of serial connection.
 - Optionally a 2nd PCV can make both reactors work at the same pressure, avoiding re-pressurization when switching reactors.






3_Parallel configuration (includes serial mode)

- Complete parallel reactor system.
- 2 MFC per reactor.
- Independent feed each reactor.
- Independent working conditions, high pressure, high temperature applications.
- 1 x HP L/G separator (L1) each reactor (total 2).
- Possibility of serial connection.
- Automatic bypass valve for L/G separator bypassing, sending all products from reactor 1 to to reactor 2.
- Optionally an automatic switching valve in the gas outlets for GC/vent selection.





OPTIONAL EQUIPMENT




Mass flow meter in the gas outlet for mass balance calculations.



Automatic six port valves for different purposes L/G separator bypass Up/down flow selection.





Fluidized bed reactor with differential pressure measurement.





High pressure 200 bar operation as an option.

mercustivity offi

8



Liquid feed


Up to 2 HPLC pumps 0.02 - 5 ml/min. Heated option up to 90 °C for heavy liquids.

Also Syringe pump for very low and accurate flow (from 0,01 μ L/min).






Up to 2 Extra **liquid evaporators temperature controlled** (50 to 450 °C) for heavier products or higher flow rates.





Gas feed

Up to 6 MFCs, 3 as standard. Easy calibration for other gases.







COMPARATIVE TABLE

	MA EFFI SOLO		MA EFFI	
	Standard	Optional	Standard	Optional
Number of reactors	1	-	1	2
HP L/G separator (L1)	-	Optional	Included	2
HP L/L/G separator (L2, i.e. FT)	-	-	-	1
Temperature	Max. 800 °C with SS316 reactor	Max. 1100 °C with special alloys	Max. 800 °C with SS316 reactor	Max. 1100 °C with special alloys
Pressure	Atmospheric	100 bar	100 bar	200 bar
Gas feed	2 MFC	Up to 4 MFC	3 MFC	Up to 6 MFC
Liquid feed	-	1 HPLC pump	-	Up to 2 HPLC pump (or syringe)
Heated liquid feed and head pump	-	1	-	2
Scale	-	1	-	2
Mass flow Meter in the gas outlet	-	1	-	2
Reactor bypass valve	Manual	Automatic	Automatic	N/A
L/G separator bypass valve	-	-	-	Optional
Up/Down flow selection valve	-	-	-	Optional
Reactor size (metallic)	9.1mm ID	Smaller (7.9mm and 5.1mm ID)	9.1mm ID	Smaller and bigger (5.1, 7.9, 13.1, 17.5 and 23.8mm ID)
Quartz reactor	-	Optional	-	Optional
Special alloys	-	Optional	-	Optional
Wax trap	-	-	-	2
Fluidized bed reactor	-	-	-	Optional
Liquid Multisampler	-	-	-	2



C/ Francisco Gervás 11

Pol. Ind. Alcobendas 28108 Alcobendas (Madrid) - SPAIN

Tel.: +34 914 840 183 - Fax: +34 916 620 141

pid@micromeritics.com

www.pidengtech.com

