

**Peer Reviewed Papers Published in 2007
Citing the ASAP 2420**

The ASAP 24** Series of Multi-Port Physisorption Analyzers



The fully automated [ASAP® 24** Accelerated Surface Area and Porosimetry Systems](#) are designed for high-performance/high sample throughput to help today's busy laboratories expand their workflow while providing highly accurate and precise surface area and porosimetry measurements. With six independently operated analysis ports, a new analysis can begin as soon as another is finished. This provides an important advantage over multi-port instruments that require all samples to be prepared or analyzed at the same time. In addition, a key, standard feature of the ASAP 24** series is a programmable sample preparation module with twelve independently operated ports. Samples may be added or removed

from degas ports without disturbing the treatment of other samples undergoing preparation or analysis. The current model is the ASAP 2420, which is available in configurations for Nitrogen, Krypton (low surface area) and Microporosity analyses.

Title	Author(s)	Publication	MIC Product(s) Cited
A bifunctional catalyst for the single-stage water-gas shift reaction in fuel cell applications. ...	KG Azzam, IV Babich, K Seshan, L Lefferts	<i>Journal of Catalysis</i> , 2007	ASAP 24** + AutoChem 29** + Pulse Chemisorb 27**
A comparative study on the dehydration of monoethanolamine over cesium phosphate modified zeolite catalysts	Sun, Y. / Yan, H. / Liu, D. / Zhao, D.	<i>Catalysis Communications</i> , 9 (5), p.924-930, Mar 2008	ASAP 24** + AutoChem 29**
A method to obtain a high loading of nano-sized Pt particles on carbon supports with a low surface area	Lisitsyn, A.S. / Okhlopkova, L.B. / Shalagina, A.E. / Ismagilov, Z.R.	<i>Carbon</i> , 46 (3), p.549-552, Mar 2008 <i>Journal of Analytical and Applied Pyrolysis</i> , 81 (1), p.33-36, Jan 2008	ASAP 24** + ASAP 2010
A study of characteristics of activated carbon produced from Jordanian olive cake	Aljundi, I.H. / Jarrah, N.		ASAP 24**
A surfactant-assisted hydrothermal deposition method for preparing highly dispersed W/?-Al₂O₃ hydrodenitrogenation...	Fan, Y., Bao, X., Wang, H., Chen, C., Shi, G.	<i>Journal of Catalysis</i> , Jan 2007 <i>Colloids and Surfaces A: Physicochemical and Engineering Aspects</i> , 316 (1), p.202-209, Mar 2008	ASAP 24**
Adsorption of anionic polyelectrolyte and comb polymers onto lead magnesium niobate	Sakar-Deliormanli, A. / Celik, E. / Polat, M.		ASAP 24** + AccuPyc 1330
Adsorption of chlorobenzene on V 2 O 5/Al 2 O 3 catalyst under dynamic conditions	LD Asn in, AA Fedorov, YS Chekryshkin	<i>Russian Journal of Applied Chemistry</i> , 2007	ASAP 24**
Bifunctional catalysts for single-stage water-gas shift reaction in fuel cell applications. Part 1. ...	KG Azzam, IV Babich, K Seshan, L Lefferts	<i>Journal of Catalysis</i> , 2007	ASAP 24** + AutoChem 29** + Pulse Chemisorb 27**
Ce-silica mesoporous SBA-15-type materials for oxidative catalysis: Synthesis, characterization, and catalytic...	Timofeeva, M.N., Jhung, S.H., Hwang, Y.K., Kim, D.K., Panchenko, V.N., Melgunov, M.S., Chesalov, Yu.A., Chang, J.S.	<i>Applied Catalysis A, General</i> , Jan 2007 <i>Microporous and Mesoporous Materials</i> , 108 (1), p.13-21, Feb 2008	ASAP 24**
Characterization and catalytic performance of SAPO-11/Hβ composite molecular sieve compared with the mechanical mixture	Zhang, X. / Wang, J. / Zhong, J. / Liu, A. / Gao, J.		ASAP 24**
Comparative characterization of polymethylsiloxane hydrogel and silylated fumed silica and silica gel	Gun'ko, V.M., Turov, V.V., Zarko, V.I., Goncharuk, E.V., Gerashchenko, I.I., Turova, A.A., Mironyuk, I.F.,.....	<i>Journal of Colloid And Interface Science</i> , Apr 2007	ASAP 24**
Composites "lithium halides in silica gel pores": Methanol sorption equilibrium	LG Gordeeva, A Freni, TA Krieger, G Restuccia, YI ...	<i>Microporous and Mesoporous Materials</i> , 2007	ASAP 24**

DC conductivity, cationic exchange capacity, and specific surface area related to chemical ...	F Henn, C Durand, A Cerepi, E Brosse, JC Giuntini	<i>Journal of Colloid And Interface Science</i> , 2007	ASAP 24**
Deactivation and Reuse of Cesium-Containing Heteropolyacid for the Isomerization of THDCPD	BH Jeong, JS Han, SW Ko, JH Lee, BJ Lee	<i>J. Ind. Eng. Chem</i> , 2007 <i>Theoretical Foundations of Chemical Engineering</i> , 2007	ASAP 24**
Dehydration of methyl phenyl carbinol catalyzed by modified active aluminum oxide	AA Lamberov, KK Gil'manov, SR Egorova, SV Trifonov		ASAP 24**
Determination of the percolation properties and pore connectivity for mesoporous solids using NMR ...	EL Perkins, JP Lowe, KJ Edler, N Tanko, SP Rigby	<i>Chemical Engineering Science</i> , 2007	ASAP 24** + AutoPore 94**
Development of a base set of toxicity tests using ultrafine TiO ₂ particles as a component of nanoparticle risk...	Warheit, D.B.; Hoke, R.A.; Finlay, C.; Donner, E.M.; Reed, K.L.; Sayes, C.M.	<i>Toxicology Letters</i> , 171 (3), p.99-110, Jul 2007	ASAP 24**
Development of porosity in carbons from yeast grains by activation with alkali metal carbonates	Urabe, Y. / Ishikura, T. / Kaneko, K.	<i>Journal of Colloid And Interface Science</i> , 319 (1), p.381-383, Mar 2008	ASAP 24**
Effect of acid/base treatment to carbon blacks on preparation of carbon-supported platinum nanoclusters	Kim, S., Park, S.J.	<i>Electrochimica Acta</i> , Feb 2007	ASAP 24**
Effect of TiOSO ₄ hydrothermal hydrolysis conditions on TiO ₂ morphology and gas-phase oxidative ...	DV Bavykin, VP Dubovitskaya, AV Vorontsov, VN ...	<i>Research on Chemical Intermediates</i> , 2007	ASAP 24**
Effects of water vapor pretreatment time and reaction temperature on CO ₂ capture characteristics of a sodium-based...	Seo, Y.; Jo, S.H.; Ryu, C.K.; Yi, C.K.	<i>Chemosphere</i> , 69 (5), p.712-718, Oct 2007	ASAP 24**
Full phenol peroxide oxidation over Fe - MMM-2 catalysts with enhanced hydrothermal stability	MN Timofeeva, MS Mel'gunov, OA Kholdeeva, ME ...	<i>Applied Catalysis B, Environmental</i> , 2007	ASAP 24**
Heterogeneous selective oxidation of formaldehyde to formic acid on V/Ti oxide catalysts: The role ...	GY Popova, TV Andrushkevich, EV Semionova, YA ...	<i>Journal of Molecular Catalysis. A, Chemical</i> , 2007	ASAP 24**
Hydroisomerization activity and selectivity of n-dodecane over modified Pt/ZSM-22 catalysts	Wang, G. / Liu, Q. / Su, W. / Li, X. / Jiang, Z. / Fang, X. / Han, C. / Li, C.	<i>Applied Catalysis A, General</i> , 335 (1), p.20-27, Feb 2008	ASAP 24**
Impact of pore structural parameters on column performance and resolution of reversed-phase ...	R Skudas, BA Grimes, E Machtejevas, V Kudirkaite, ...	<i>Journal of Chromatography A</i> , 2007	ASAP 24**
In situ ATR-IR study of CO adsorption and oxidation over Pt/Al ₂ O ₃ in gas and aqueous phase: Promotion effects by water...	Ebbesen, S.D., Mojet, B.L., Lefferts, L.	<i>Journal of Catalysis</i> , Feb 2007	ASAP 24**
In-situ synthesis of nanoparticles via supersolubilizing micelle self-assembly	DC Wang	<i>Science in China Series B: Chemistry</i> , 2007	ASAP 24**
Large pore heavy oil processing catalysts prepared using colloidal particles as templates	Chen, S.L.; Dong, P.; Xu, K.; Qi, Y.; Wang, D.	<i>Catalysis Today</i> , 125 (3), p.143-148, Jul 2007	ASAP 24**
Liquid chromatography method for quantification of surface connected mesoporosity in ultrastable Y ...	L Teyssier, M Thomas, C Bouchy, JA Martens, E ...	<i>Microporous and Mesoporous Materials</i> , 2007	ASAP 24** + AutoPore 95**
Low-Temperature Adsorption/Storage of Hydrogen on FAU, MFI, and MOR Zeolites with Various Si/Al ...	SH Jung, JW Yoon, JS Lee, JS Chang	<i>CHEMISTRY-WEINHEIM EUROPEAN JOURNAL</i> , 2007	ASAP 24** + TriStar 3000
Macroporous molecularly imprinted polymer/cryogel composite systems for the removal of endocrine disrupting trace...	Le Noir, M.; Plieva, F.; Hey, T.; Guieysse, B.; Mattiasson, B.	<i>Journal of Chromatography A</i> , 1154 (1), p.158-164, Jun 2007	ASAP 24**
Mesoporous micelle templated silica with incorporated C 8 and C 18 phase	J Goworek, W Stefaniak, A Kierys, M Iwan	<i>Journal of Thermal Analysis and Calorimetry</i> , 2007	ASAP 24**
Methanol conversion on SAPO-34 catalysts prepared by mixed template method	Lee, Y.J.; Baek, S.C.; Jun, K.W.	<i>Applied Catalysis A, General</i> , 329, p.130-136, Oct 2007	ASAP 24**
Phase-selective crystallization of cobalt-incorporated aluminophosphate molecular sieves with large pore by microwave...	Jhung, S.H. / Jin, T. / Kim, Y.H. / Chang, J.S.	<i>Microporous and Mesoporous Materials</i> , 109 (1), p.58-65, Mar 2008	ASAP 24**
Photocatalytic decomposition of orange II over TiO ₂ -loaded on SBA-15 prepared using a microwave ...	WY Jung, BH Noh, SH Baek, GD Lee, SS Park, SS Hong	<i>Reaction Kinetics and Catalysis Letters</i> , 2007	ASAP 24**

Pore curvature and support composition effects on the electronic properties of supported Pt ...	M Lepage, T Visser, AMJ van der Eerden, F ...	<i>Vibrational Spectroscopy</i> , 2007	ASAP 24**
Production of highly mesostructured SBA-15 silicas at pH around the PZC	D Zhao, S Qiu, Y Tang, C Yu	<i>Recent Progress in Mesostructured Materials: Proceedings of ...</i> , 2007	ASAP 24**
Propene epoxidation over Au/Ti-SBA-15 catalysts	Sacaliuc, E.; Beale, A.M.; Weckhuysen, B.M.; Nijhuis, T.A.	<i>Journal of Catalysis</i> , 248 (2), p.235-248, Jun 2007	ASAP 24**
Pulmonary toxicity study in rats with three forms of ultrafine-TiO₂ particles: Differential responses related to...	Warheit, D.B., Webb, T.R., Reed, K.L., Frerichs, S., Sayes, C.M.	<i>Toxicology</i> , Jan 2007	ASAP 24**
Realumination of dealuminated HZSM-5 zeolite by citric acid treatment and its application in preparing FCC gasoline...	Fan, Y., Lin, X., Shi, G., Liu, H., Bao, X.	<i>Microporous and Mesoporous Materials</i> , Jan 2007	ASAP 24** ASAP 24** +
Single stage water gas shift conversion over Pt/TiO₂—Problem of catalyst deactivation	KG Azzam, IV Babich, K Seshan, L Lefferts Ramachandran, C.E. / Du, H. / Kim, Y.J. / Kung, M.C. / Snurr, R.Q. / Broadbelt, L.J.	<i>Applied Catalysis A, General</i> , 2007	Pulse Chemisorb 27**
Solvent effects in the epoxidation reaction of 1-hexene with titanium silicalite-1 catalyst		<i>Journal of Catalysis</i> , 253 (1), p.148-158, Jan 2008	ASAP 24**
Studies of Iron Effects on FCC Catalysts	S CangZhou Refinery, C HeBei	<i>Fluid Catalytic Cracking VII: Materials, Methods and Process ...</i> , 2007	ASAP 24** + AutoChem 29**
ŠTÚDIUM POVRCHOVÝCH VLASTNOSTÍ SLOVENSKÉHO PRÍRODNÉHO ZEOLITU–KLINOPTILOLITU FYZIKÁLNOU ADSORPCIOU ...	F Mária, H Pavol	<i>Petroleum & Coal</i> , 2007	ASAP 24**
Study of palladium catalyst deactivation in synthesis of 4, 10-diformyl-2, 6, 8, 12-tetraacetyl-2, 4 ...	AP Koskin, IL Simakova, VN Parmon	<i>Reaction Kinetics and Catalysis Letters</i> , 2007	ASAP 24**
Surface characterisation of chemically reduced electrolytic manganese dioxide	Malloy, A.P. / Donne, S.W. Gun'ko, V.M.; Blitz, J.P.; Gude, K.; Zarko, V.I.; Goncharuk, E.V.; Nychiporuk, Y.M.; Lebeda, R.; (...); Ptushinskii, Y.G.	<i>Journal of Colloid And Interface Science</i> , 320 (1), p.210-218, Apr 2008	ASAP 24**
Surface structure and properties of mixed fumed oxides		<i>Journal of Colloid And Interface Science</i> , 314 (1), p.119-130, Oct 2007	ASAP 24**
Synergistic effect of polymer-surfactant mixtures on the stability of aqueous silica suspensions	A Sakar-Deliormanli	<i>Journal of the European Ceramic Society</i> , 2007	ASAP 24**
Synthesis of aluminum oxides from the products of the rapid thermal decomposition of hydrargillite ...	YY Tanashev, EM Moroz, LA Isupova, AS Ivanova, GS ...	<i>Kinetics and Catalysis</i> , 2007	ASAP 24**
Synthesis of LiCo_{0.3}Ni_{0.7}O₂ as cathode materials for lithium ion batteries by oxidation-ion exchange of...	Tong, D.G. / Tang, A.D. / Chu, W. / Tang, L.X. / Huang, K.L. / He, Y. / Ji, X.Y.	<i>Materials Chemistry & Physics</i> , 107 (2), p.385-391, Feb 2008	ASAP 24**
Synthesis of nano-sized porous γ-alumina powder via a precipitation/digestion route	HS Potdar, KW Jun, JW Bae, SM Kim, YJ Lee	<i>Applied Catalysis A, General</i> , 2007	ASAP 24**
Synthesis of Nanosize-Necked Structure α- and γ-Fe₂O₃ and its Photocatalytic Activity	SK Apte, SD Naik, RS Sonawane, BB Kale, JO Baeg	<i>Journal of the American Ceramic Society</i> , 2007	ASAP 24**
Synthesis of Nanosize-Necked Structure-α-Fe₂O₃ and its Photocatalytic Activity	SK Apte, SD Naik, RS Sonawane, BB Kale, JO Baeg	<i>Journal of the American Ceramic Society</i> , 2007	ASAP 24**
Synthesis of thermo-stable high surface area alumina powder from sol-gel derived boehmite	Kim, S.M.; Lee, Y.J.; Jun, K.W.; Park, J.Y.; Potdar, H.S.	<i>Materials Chemistry & Physics</i> , 104 (1), p.56-61, Jul 2007	ASAP 24**
Temperature Study of Nonionic Polymers Adsorption at the Alumina-Solution Interface	M Wisniewska	<i>Journal of the American Ceramic Society</i> , 2007	ASAP 24**
Textural properties vs. CEC and EGME retention of Na-X zeolite prepared from fly ash at room temperature	Derkowski, A., Franus, W., Waniak-Nowicka, H., Czimerova, A.	<i>International Journal of Mineral Processing</i> , Mar 2007	ASAP 24**
The modification and stability of γ-Al₂O₃ based catalysts for hydrolytic decomposition of CF₄	Xu, X.F., Jeon, J.Y., Choi, M.H., Kim, H.Y., Choi, W.C., Park, Y.K.	<i>Journal of Molecular Catalysis. A, Chemical</i> , Apr 2007	ASAP 24**
The synthesis of mesoporous aluminosilicate using microcline for adsorption of mercury(II)	Wu, X.W. / Ma, H.W. / Li, J.H. / Zhang, J. / Li, Z.H.	<i>Journal of Colloid And Interface Science</i> , 315 (2), p.555-561, Nov 2007	ASAP 24**

The temperature influence on the adsorption and electrokinetical properties in the nonionic polymer/ ...	M Wisniewska, S Chibowski, T Urban	<i>Materials Chemistry & Physics, 2007</i>	ASAP 24**
Toward commercialization of nanozeolites	R Xu, Z Gao, J Chen, W Yan	<i>From Zeolites to Porous MOF Materials-the 40th Anniversary ..., 2007</i>	ASAP 24**
Two-Dimensional Nanosheet Crystals* Water and Solute Movement in Air Pollution Control Residues Measured by Time-Domain Reflectometry	J Seo, Y Jun, S Park, H Nah, T Moon, B Park, JG ...	<i>Angew. Chem. Int. Ed, 2007</i>	ASAP 24**
Y/composite titania-silica (CTS) supported catalyst for hydrotreating coker gas oil	M van Praagh, M Persson, KM Persson	<i>Journal of Environmental Engineering, 2007</i>	ASAP 24**
ZrFe INTERMETALLIDES FOR FISCHER-TROPSCH SYNTHESIS: PURE AND ENCAPSULATED INTO ALUMINA-CONTAINING ...	Zhou, Y.; Wei, Q.; Ma, H.; Zhang, Z.	<i>Catalysis Today, 125 (3), p.211-219, Jul 2007</i>	ASAP 24**
???? ? ? ? ? ? ? ? ?	BH Davis, ML Occelli	<i>Fischer-Tropsch Synthesis, Catalysts and Catalysis, 2007</i>	ASAP 24**
???? ? ? ? ? ? ? ? ?	? ? ? , ? ? ? , ? ? ? ,	? ? ? ? ? ? ? ? , 2007	ASAP 24**
???? ? ? ? ? ? ? ? ?	? ? ? , ? ? ? ?	... Model X'Pert-MPD,	
???? ? ? ? ? ? ? ? ?	? ? ? , ? ? ? ?	<i>Phillips Analytical</i>	
???? ? ? ? ? ? ? ? ?	? ? ? ? ? ? ? ? ? ? 2007?	? ? ? ? FS?	
???? ? ? ? ? ? ? ? ?	? ? ? ? ? ? ? ? ? ? ? ? ? ,	? ? ? ? ? ? ? ? ? ?	ASAP 24**
Tetrahydrodicyclopentadiene?	? ? ? ?	<i>J. Korean Ind. Eng. Chem, 2007</i>	ASAP 24**
? ? ? ? ? ? ? ? ? ? ? ? ? ?	? ? ? ?	? ? ? ? ? ? ? ? ? ? ,	
? ? ? ?	? ? ? , ? ? ? ? , ? ? ? ?	2007 ? ? ? ? ?	ASAP 24**
? ? ? ? NiB/MCM-41? ? ? ? 3-(N-?	? ? ? , ? ? ? ? , ? ? ? ? ,		
? ? ? ? 4? ? ? ? ? ? ? ?	? ? ? , ? ? ? ? , ...	? ? ? ? , 2007 ? ? ? ? ?	ASAP 24**

An additional 119 papers citing the ASAP 24** series can be found in ***Bibliography of Peer-reviewed Papers Citing Micromeritics Products***, a downloadable file located in the Applications section of Micromeritics' web site (http://www.micromeritics.com/applications/bibliography_of_papers.aspx)