

5. Ordering Information

The FlowPrep 060 System components and accessories can be ordered using one of the following methods:

- Call our Customer Service Department at (770) 662-3636
- Access our web site at www.micromeritics.com
- Contact your local sales representative

When ordering, please use the information provided below to place your order

Table 5-1. FlowPrep Components and Accessories

| Part Number | Item and Description |
|--------------|---|
| 236-33005-00 | PrepSeal Kit, used to degas and transfer air-sensitive samples to an analysis port without contamination. |
| 061-33001-00 | Flow Control Valve Kit; allows you to regulate the flow rate to the sample tube |
| 060-25809-00 | Gas Delivery Tube, stainless steel, set of 6 |
| 060-25804-00 | Flexible Tubing, internal (valve inlet), set of 6 |
| 060-25805-00 | Flexible Tubing, external, set of 6 |
| 003-51147-00 | Slow-Blow Fuse, 2.5 A, 100-200 VAC |
| 003-51148-00 | Slow-Blow Fuse, 1.25 A, 220-240 VAC |
| 003-51197-00 | Slow-Blow Fuse, 1.25 A, 220-240 VAC (Type T, time lag) |
| 230-02001-00 | Gas inlet line, 1/8 in. x 6 ft. long, copper |
| 004-32004-00 | Sample Tube Stopper, rubber |
| 004-62230-58 | Gas Regulator, for use with helium or nitrogen, 30 psig |
| 004-25549-00 | Reducer Fitting, 1/4-in. for 1/8-in. copper tubing |
| 065-25815-00 | Heating Station Adapter, TriStar and ASAP Series analyzers |
| 065-25811-00 | Heating Station Adapter, Gemini bulb sample tube, 3/4 in |
| 065-25812-00 | Heating Station Adapter, Gemini straight-wall sample tube, 3/8 in. |
| 003-22006-01 | Temperature Controller |
| 060-42803-00 | Operator's Manual |