



GEOPYC 1360

SPECIFICATIONS

Reproducibility:	When sample volume is at least 25% of sample holder volume: Typically $\pm 1.1\%$
Sample Chambers:	<p>12.7-mm ID: sample volumes between 0.3 and 0.8 cm³</p> <p>19.1-mm ID: sample volumes between 0.8 and 2.4 cm³</p> <p>25.4-mm ID: sample volumes between 2.4 and 5.3 cm³</p> <p>38.1-mm ID: sample volumes between 5.3 and 13 cm³</p> <p>50.8-mm ID: sample volumes between 13 and 25 cm³</p>
Physical:	<p>Width: 22 in</p> <p>Height: 11 in</p> <p>Depth: 15 in</p> <p>Weight: 42 lbs</p>
Electrical Power:	95 VA, Voltage: 85 to 265 VAC, Frequency: 47 to 63 Hz
Miscellaneous:	<p>ISO 9001 manufacturer</p> <p>CE certified</p>