

SPECIFICATIONS

AccuPyc[®] 1340 - TEC

Thermo Electric Temperature Control



TEMPERATURE CONTROL

Range	15 °C to 36 °C, (+/- 0.1 °C)
Control	Adjustable in 0.1 °C increments
Environment	Stable between 15 °C to 35 °C (59 to °F) 20% to 80% relative humidity

PHYSICAL

Height	27.3 cm (10.7 in)
Width	28.7 cm (11.3 in)
Depth	36.1 cm (14.2 in.)
Weight	13.8 kg (30.5 lbs)

ELECTRICAL

Voltage	90 to 264 VAC
Power	100 VA MAX
Frequency	50 to 60 Hz

SAMPLE CHAMBER

10 cm³ chamber Available disposable insert cups	1.80 cm ID x 3.93 cm D (0.72 in. ID x 1.55 in. D.)
--	--

Due to continuous improvements, specifications are subject to change without notice.