

Introducing...




(patent pending)


The innovative *SPECTester* measures your sample, of **up to 6 unique components** and, with a touch of the finger, reports **how much** as well as **why** your material is segregating. Fully automated, it provides data about: component concentrations, particle size differences, product uniformity, and **up to 4 specific segregation mechanisms**. This segregation tester can be used in R&D facilities as well as in production plants for on-the-spot, mid-stream quality control.

Materials of different sizes, shapes, colors and/or chemical composition

The SPECTester handles them all



Touch Screen Data Entry
Fast Graphical and Numerical Results
UBS Thumb-drive (included)
for Reporting Excel-ready data



3 component Mixture



THE
*SPECTester*TM

