


Sub-Sieve Size and Area Report

Logo not found

Instrument: MIC SAS	Serial number: 5397	Software: MIC SAS II 5800 v1.00
Started: Jan. 27, 2021, 6:42 p.m.	Ended: Jan. 27, 2021, 6:51 p.m.	
Material: Portland Cement	Reported: Feb. 4, 2021, 10:55 a.m.	
Sample name: Portland Cement	Operator: JGS	
Batch name: 46h NIST	Collate results: No	
Mass: 3.1490 g	Absolute density: 3.1500 g/cm ³	



Porosity	Particle Size (µm)	Surface Area (cm ² /g)	Fisher Particle Size (µm)	Fisher Surface Area (cm ² /g)
0.424	4.821	3951	4.571	4167